

*Angelman Syndrome Foundation Communication Training Series  
December 8, 2015*

## Getting Started with Proloquo2Go

Jennifer Marden, MA CCC-SLP  
AssistiveWare

## Agenda

- ❖ Introduction
- ❖ Creating a new Crescendo Intermediate Core user
- ❖ Getting acquainted with your new user
- ❖ Personalizing your new user
- ❖ Customizing for physical and visual access
- ❖ Customizing for literacy

## For More Information

- ❖ [www.assistiveware.com/e-learning-videos](http://www.assistiveware.com/e-learning-videos)
  - ❖ Getting started with Proloquo2Go 4
  - ❖ Transitioning from Proloquo2Go 3 to Proloquo2Go 4
  - ❖ Basic Customization
  - ❖ Literacy Development
  - ❖ Customizing for vision
  - ❖ Access Methods - Hold Duration and Select on Release
- ❖ [bit.ly/Proloquo2Go-Family](http://bit.ly/Proloquo2Go-Family)
- ❖ [bit.ly/Proloquo2Go-Professionals](http://bit.ly/Proloquo2Go-Professionals)
- ❖ [bit.ly/howfindwords](http://bit.ly/howfindwords)

## Moving to Crescendo Core

- ❖ Back up your current user!
- ❖ Update to Proloquo2Go 4.1.1, iOS 8.2 at least
- ❖ Create a new Crescendo Intermediate Core user, 6x10 grid size at least
- ❖ Copy any customizations you need from original user
- ❖ Personalize vocabulary
- ❖ Back up your new user!